

**1. Part Number:** A335-012

**2. Photograph:**



### 3. Product Information:

Description	Volumetric Spoon 1.25ml
Material of Construction:	304 stainless steel
Method of Construction:	Single piece pressing
Nominal Volume:	1.25ml
Nominal Weight:	11g
Bowl Diameter:	23mm
Overall Length:	108mm
Surface Finish:	Better than 0.5 microns Ra
Recommended Storage Conditions:	Dry and ambient temperature

Note: The volumes of these measuring spoons are not certified.

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24