

1. Part Number: A335-050**2. Photograph:****3. Product Information:**

Description	Volumetric Spoon 5ml
Material of Construction:	304 stainless steel
Method of Construction:	Single piece pressing
Nominal Volume:	5ml
Nominal Weight:	16g
Bowl Diameter:	33mm
Overall Length:	112mm
Surface Finish:	Better than 0.5 microns Ra
Recommended Storage Conditions:	Dry and ambient temperature

Note: The volumes of these measuring spoons are not certified.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24