

1. Part Number: A335-1000

2. Photograph:



3. Product Information:

| | |
|---------------------------------|-----------------------------|
| Description | Volumetric Spoon 100ml |
| Material of Construction: | 304 stainless steel |
| Method of Construction: | Single piece pressing |
| Nominal Volume: | 100ml |
| Nominal Weight: | 52g |
| Bowl Diameter: | 79mm |
| Overall Length: | 127mm |
| Surface Finish: | Better than 0.5 microns Ra |
| Recommended Storage Conditions: | Dry and ambient temperature |

Note: The volumes of these measuring spoons are not certified.

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 18-Dec-24