

**1. Part Number:** A335-150

**2. Photograph:**



### 3. Product Information:

|                                 |                             |
|---------------------------------|-----------------------------|
| Description                     | Volumetric Spoon 15ml       |
| Material of Construction:       | 304 stainless steel         |
| Method of Construction:         | Single piece pressing       |
| Nominal Volume:                 | 15ml                        |
| Nominal Weight:                 | 25g                         |
| Bowl Diameter:                  | 47mm                        |
| Overall Length:                 | 119mm                       |
| Surface Finish:                 | Better than 0.5 microns Ra  |
| Recommended Storage Conditions: | Dry and ambient temperature |

Note: The volumes of these measuring spoons are not certified.

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24