

1. Part Number: A335-1

2. Photograph:



3. Product Information:

Description Stainless Measuring Spoons (1 set of 4 spoons)
Material of Construction: 304 stainless steel
Method of Construction: Single piece pressing
Nominal Volume: The set comprises of 4 spoons:

Volume	Overall Length	Max. Bowl Diameter	Weight
1.25ml	107mm	23mm	12g
5ml	112mm	33mm	17g
15ml	118mm	46mm	25g
100ml	127mm	79mm	55g

Note: The volumes of these measuring spoons are not certified.

Surface Finish: Better than 0.5 microns Ra
Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 18-Dec-24