

1. Part Number: A335-200

2. Photograph:



3. Product Information:

Description	Volumetric Spoon 20ml
Material of Construction:	304 stainless steel
Method of Construction:	Single piece pressing
Nominal Volume:	20ml
Nominal Weight:	31g
Bowl Diameter:	51mm
Overall Length:	141mm
Surface Finish:	Better than 0.5 microns Ra
Recommended Storage Conditions:	Dry and ambient temperature

Note: The volumes of these measuring spoons are not certified.

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 18-Dec-24