

1. Part Number: A346-1500

2. Photograph:



3. Product Information:

Description	Scorpion PharmaScoop® 1500ml
Nominal Volume:	1500ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Body Length:	230mm
Body Width:	155mm
Maximum Height (inc. handle):	140mm
Weight:	630g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Sep-21