


Specification Sheet

Specification Number: A418-1
Revision: 03
Date Issued: 27-Oct-20

1. Part Number: A418-1

2. Photograph:





3. Product Information:

Description Stainless Container 1L with Internal GL45 Thread

Materials of Construction:
- Body: 316L stainless steel
- Lid: 316L stainless steel

- Lid: 316L stainless stee - Clamp: 304 stainless steel

- Gasket: Silicone (FDA acceptable grade)

- Cap: 316L stainless steel

- Cap Seal: Silicone (FDA acceptable grade)

Method of Construction:

All welds ground and polished

Crevice free interior

Surface Finish: Better than 0.5 microns Ra

Nominal Volume 1 Litre

Dimensions A = 119mm

B = 101mm C = 152mm D = 104mm

Weight: 1.6 Kg (body + lid + gasket + clamp)

Recommended Storage Conditions: Dry and ambient temperature

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 27-Oct-20



4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com