


1. Part Number: A454-50

2. Photograph:



3. Product Information:

Description:	Storage Container 50L
Material of Construction:	304 Stainless Steel (body and lid)
Method of Construction:	Side Handles are spot welded to body Patch plate for lid handle is spot welded to lid
Nominal Volume:	50L
Dimensions:	A = 560 mm B = 535 mm C = 435 mm D = 350 mm E = 355 mm
Graduation Marks:	Every 10L
Handle Type:	Side handles
Nominal Weight:	7.45 kg (body & lid)
Surface Finish:	Bright Polish (Ra < 0.5 microns)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 6-Jan-25