

**1. Part Number:** A516-120

**2. Photograph:**



### 3. Product Information:

Description	Sieving Tool 120mm
Material of Construction:	316L stainless steel
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Base Diameter:	120mm
Base Thickness:	10mm
Overall Height:	160mm
Nominal Weight:	1130g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Storage Conditions:	Dry, ambient conditions

### 4. Approval:

Signed:



**Name:** Kuldeep Flora  
**Position:** Quality Manager  
**Date:** 6-Jan-25