

1. Part Number: A534-2000

2. Photograph:



3. Product Information:

Description	Stainless Mixing Bowl 2000ml
Nominal Volume:	2000ml
Material of Construction:	304 stainless steel
Method of Construction:	Single piece construction
Surface Finish:	Better than 0.5 microns Ra
Height	112mm
Top Diameter:	217mm
Weight:	238g
Graduation marks:	None
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 6-Jan-25