

**1. Part Number:** A618-235

**2. Photograph:**



### 3. Product Information:

Description:	Narrow Spoon Spatula 235mm
Material of Construction:	304 Stainless steel
Overall Length:	235mm
Bowl Length:	40mm
Bowl Width:	18mm
Nominal Bowl Volume:	2ml
Nominal Weight:	22g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 7-Jan-25