


1. Part Number: A622-310**2. Photograph:****3. Product Information:**

Description:	Non-Slip Beaker Tongs 310mm
Material of Construction:	- 304 Stainless steel - PVC coating on end of arms
Overall Length:	310mm
Diameter of Arms:	6mm
Max. Height of Arms:	22mm
Nominal Weight:	219g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

Note: Suitable for use with beakers over 45mm in diameter

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Jan-25