


1. Part Number: A626-10

2. Photograph:



A = 295 mm
 B = 280 mm
 C = 255 mm
 D = 220 mm
 E = 225 mm

3. Product Information:

Description:	Mini-Drum with Toggles - 10L
Method of Construction:	- Body welds ground and polished - Patch plates for toggles, swing handle and lid handle are spot welded
Materials of Construction:	- Body: 316L stainless steel - Lid: 316L stainless steel - Toggle Clamp: 304 stainless steel - O-ring: Silicone Rubber (conforms to FDA 177.2600)
Nominal Volume:	10L
Surface Finish:	Bright Polish (Ra < 0.5 microns)
Graduation Marks:	None
Handle Type:	Swing Handle
Number of Toggles	3
Nominal Weight:	2.6 kg (body + lid + toggles + O-ring)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 9-Jul-21