


Specification Sheet

Specification Number: A626-12 Revision: 03 Date Issued: 9-Jul-21

1. Part Number: A626-12

2. Photograph:





3. Product Information:

Description: Mini-Drum with Toggles - 12L

Method of Construction: - Body welds ground and polished

- Patch plates for toggles, swing handle and lid handle are spot

welded

Materials of Construction: - Body: 316L stainless steel

> 316L stainless steel - Lid: - Toggle Clamp: 304 stainless steel

- O-ring: Silicone Rubber (conforms to FDA 177.2600)

Nominal Volume:

Surface Finish: Bright Polish (Ra < 0.5 microns)

Graduation Marks: None

Handle Type: Swing Handle

Number of Toggles 3

Nominal Weight: 3.0 kg (body + lid + toggles + O-ring)

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:

Kuldeep Flora Name:

Position: **Quality Manager**

Date: 9-Jul-21

