

Specification Sheet Specification Number: A626-1 Revision: 03


Date Issued:

03 9-Jul-21

- 1. Part Number:
- A626-1

2. Photograph:





## 3. Product Information:

Description: Method of Construction:

Materials of Construction:

Nominal Volume:

Graduation Marks:

Number of Toggles

**Recommended Storage Conditions:** 

Nominal Weight:

4. Approval:

Surface Finish:

Handle Type:

Mini-Drum with Toggles - 1L

- Body welds ground and polished
- Patch plates for toggles, swing handle and lid handle are spot welded
- Body: 316L stainless steel
- Lid: 316L stainless steel
- Toggle Clamp: 304 stainless steel
- O-ring: Silicone Rubber (conforms to FDA 177.2600)
- 1L

Bright Polish (Ra < 0.5 microns)

None

Swing Handle

3

0.65 kg (body + lid + toggles + O-ring)

Dry and ambient temperature

Signed:

Name:Kuldeep FloraPosition:Quality ManagerDate:9-Jul-21



4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com