


1. Part Number: A626-6

2. Photograph:



A = 280 mm
B = 270 mm
C = 205 mm
D = 175 mm
E = 180 mm

3. Product Information:

Description: Mini-Drum with Toggles - 6L

Method of Construction:

- Body welds ground and polished
- Patch plates for toggles, swing handle and lid handle are spot welded

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Toggle Clamp: 304 stainless steel
- O-ring: Silicone Rubber (conforms to FDA 177.2600)

Nominal Volume: 6L

Surface Finish: Bright Polish (Ra < 0.5 microns)

Graduation Marks: None

Handle Type: Swing Handle

Number of Toggles: 3

Nominal Weight: 1.8 kg (body + lid + toggles + O-ring)

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 9-Jul-21