


1. Part Number: A645-12

2. Photograph:



3. Product Information:

Description:	Test Tube Rack (2x6)
Material of Construction:	Stainless steel
Overall Length:	152mm
Width:	58mm
Height:	90mm
Number of Holes:	12 (2x6)
Hole Diameter:	19mm
Nominal Weight:	220g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 8-Jan-25