

**1. Part Number:** A716-300

**2. Photograph:**



**3. Product Information:**

Description:	Powder Funnel 300mm
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	8 Litres
Top Diameter:	300mm
Top Diameter including flange:	330mm
Overall Height:	310mm
Spout Length:	25mm
Spout Diameter:	Spout ends with a 2" ferrule (BS4825-3)
Nominal Weight:	1720g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25