

1. Part Number: A725-300

2. Photograph:



3. Product Information:

Description:	Stainless Bucket 300ml
Material of Construction:	Stainless steel (grade of stainless steel not specified)
Surface Finish:	Better than 0.5 microns Ra, crevice free interior of bucket
Construction:	Handle support spot welded to side of bucket
Nominal Volume:	300 mL
Brimful Volume:	350 mL
Top Diameter of Bucket:	92 mm
Overall Height of Bucket:	88 mm
Max. Width (inc. Handle)	115 mm
Max. Height (inc. Handle)	135 mm
Graduation Marks:	No
Nominal Weight:	105 g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25