


1. Part Number: A836-1000

2. Photograph:



3. Product Information:

Description: Stainless Bottle 1000 ml

Material of Construction:

Body:	316L stainless steel
Lid:	316L stainless steel
Clamp:	304 stainless steel
Gasket:	Silicone (conforms to FDA CFR 177.2600)

Method of Construction: Crevice free, all interior welds ground and polished

Surface Finish: Better than 0.5 microns Ra

Nominal Volume: 1000ml

Brimful Volume: 1200ml

Dimensions:

A = 64 mm
B = 48 mm
C = 178 mm
D = 102 mm

Nominal Weight (body, lid, clamp & gasket): 980 g

Recommended Storage Conditions: Dry and ambient temperature

BSE/TSE statement: All polishing compounds used in the manufacture of this scoop are of vegetable origin

Note: The neck of the container is fitted with a 2" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Dec-24