


1. Part Number: A886-1000

2. Photograph:



3. Product Information:

Description	Stainless Bottle GL45 - 500ml
Materials of Construction:	- Body: 316L stainless steel - Lid: 316L stainless steel - Gasket: PTFE (FDA 21 CFR 177.1550 compliant)
Nominal Volume:	1000ml
Brimful Volume:	1260ml
Method of Construction:	Crevice free interior, welds ground & polished
Surface Finish:	Electropolished (internal and external) Better than 0.5 microns Ra
Dimensions:	A = 50 mm B = 28 mm C = 175 mm D = 108 mm
Nominal Weight (body, lid & gasket):	660g
BSE/TSE statement:	BSE/TSE Free (the electropolishing process removes all polishing waxes)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 14-Jan-25