


1. Part Number: A886-2500

2. Photograph:




3. Product Information:

Description	Stainless Bottle GL45 - 2500ml
Materials of Construction:	<ul style="list-style-type: none"> - Body: 316L stainless steel - Lid: 316L stainless steel - Gasket: PTFE (FDA 21 CFR 177.1550 compliant)
Nominal Volume:	2500ml
Brimful Volume:	2770ml
Method of Construction:	Upper and lower sections welded. Interior weld left as laid but electropolished.
Surface Finish:	Electropolished (internal and external) Better than 0.5 microns Ra
Dimensions:	A = 50 mm B = 28 mm C = 221 mm D = 140 mm
Nominal Weight (body, lid & gasket):	930g
BSE/TSE statement:	BSE/TSE Free (the electropolishing process removes all polishing waxes)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 15-Dec-25