


**1. Part Number:** A886-500

**2. Photograph:**



### 3. Product Information:

**Description**

Stainless Bottle GL45 - 500ml

**Materials of Construction:**

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Gasket: PTFE (FDA 21 CFR 177.1550 compliant)

**Nominal Volume:**

500ml

**Brimful Volume:**

680ml

**Method of Construction:**

Crevice free interior, welds ground & polished

**Surface Finish:**

Electropolished (internal and external)  
Better than 0.5 microns Ra

**Dimensions:**

A = 50 mm  
B = 28 mm  
C = 141 mm  
D = 88 mm

**Nominal Weight (body, lid & gasket):**

507g

**BSE/TSE statement:**

BSE/TSE Free (the electropolishing process removes all polishing waxes)

**Recommended Storage Conditions:**

Dry and ambient temperature

### 4. Approval:

**Signed:**



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 14-Jan-25