

1. Part Number: A976-100

2. Photograph:



3. Product Information:

| | |
|---------------------------------|---|
| Description | Angled Forceps 100mm |
| Material of Construction: | 316L stainless steel |
| Method of Construction: | Spot welded and polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Overall Length: | 104mm |
| Nominal Weight: | 11g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are free of animal products |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 8-Jan-25