

**1. Part Number:** A986-150

**2. Photograph:**



### 3. Product Information:

Description:	Liquid Funnel 150mm
Material of Construction:	316L Stainless Steel
Top Diameter:	150mm
Outer Diameter of Spout:	20mm
Overall Length:	171mm
Nominal Spout Length:	95mm
Nominal Weight:	167g
Side Handle:	Spot welded to cone section of funnel
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



**Name:** Kuldeep Flora  
**Position:** Quality Manager  
**Date:** 8-Jan-25